ALCOVEST / ALCOVEST RAPID

Descripción

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ALCOVEST / ALCOVEST RAPID

Description

Special polysilicate and alcohol base investment for Chrome Cobalt casting on moulds.

Instructions for use

WORKING CONDITIONS

Recommended powder and liquid working temperature: 21 °C. Higher temperatures both room and of the mixing liquid accelerate the setting procedure.
Do not use crystallized liquid!!

PREPARING THE LIQUID

Add the contents of the small bottle (B) to the large one (A).
Shake it for 1 minute and leave the mixture to rest for 4 hours, after which time will be ready to use.
The mixed liquids will begin to react after 10 minutes. An increase of the temperature must be observed.
Storage: If you do not use all the liquid you have prepared within 10 days, we recommend storing it in the refrigerator, otherwise it might crystallize.

MIXING

The powder and liquid should be mixed by putting the liquid in a suitable container first before adding the powder and pasting it on the vibrator once it has been turned on.

Mixture proportion: 15 - 18 ml liquid / 100 g powder

Mixing time: 60 sec (by hand)

WORKING TIME

8 - 10 minutes / 6 minutes for RAPID version (20 - 21° C)
The working time can be reduced when the liquid is at a higher temperature, at temperatures over 21 °C, and when using contaminated utensils. The working time can be extended when the liquid and/or powder are at temperatures below 20° C.

PREPARING THE DUPLICATES

Clean the duplicate beforehand and place it on the vibrator, apply the investment with a gentle vibration and keep it on for 10 minutes. After 30-40 minutes it can be withdrawn from the gelatin

FILLING THE CYLINDERS

To fill the cylinder, prepare a slightly more liquid investment than for the duplicates, to make the filling operation easier.
Once the cylinder is dry, cut off the top layer where the superfluous liquid had collected, to prevent the cylinder from cracking when it is heated.
Setting time: 30-45 minutes.
Note: Before filling the cylinder, paint the wax mould with a fine investment, or apply a wax degreasing varnish

HEATING

After drying and cutting, place the cylinder in the cold furnace, dry it at 300°C for 30-60 minutes depending on the size of the cylinder.

Increase the temperature by between 7 and 9 degrees per minute until it reaches:

950 ° C for partial plates.

1050°C for complete plates.

Maintain the final temperature 30 – 60 minutes depending on the size and number of cylinders.

After setting, let the cylinder cool at room temperature.
Do not immerse in cold water!!

PRECAUTION

Do not mix the ALCOVEST investment in a container in which gypsum products are pasted.

NOTE: The proportion in the powder+liquid kit is thought to an easier use for the consumer. In the practical use in the laboratory, it's possible that you need more powder or liquid. To solve this, replacement for both liquid and powder are sold (see below: packaging)

Storage and expiry

The stability of the product is 4 years (3 years for the liquid) respecting the storage instructions.
Keep the packages tightly closed in a dry place.

Packaging

6630-200 ALCOVEST 20 Kg (Powder + Liquid A & B)
6630-201 ALCOVEST RAPID 20 Kg (Powder + Liquid A & B)

6630-020 ALCOVEST Replacement Liquid A & B
6630-021 ALCOVEST Replacement Liquid A & B RAPID
6630-220 ALCOVEST 20 Kg Powder

Safety information

Avoid the formation and breathing of the dust. (It contains silica)
Work in well-ventilated areas or with aspiration

Guarantee:

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet of ALCOVEST that can be requested from your usual supplier.

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