# VarseoSmile Crown <sup>plus</sup>

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**en** Instructions for use

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
Temperaturbegrenzung

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## VarseoSmile Crown<sup>plus</sup>

Resin for 3D printing of permanent single crowns, inlays, onlays and veneers.

### 1. Intended use/Indication

VarseoSmile Crown<sup>plus</sup> is a light-curing, free-flowing plastic based on methacrylic acid esters for the production of permanent single crowns, inlays, onlays and veneers.

### 2. Contraindications

Known allergy to one or more ingredients. In cases of doubt, the allergy should be clarified and ruled out based on a specific test prior to the application of this product. VarseoSmile Crown<sup>plus</sup> should not be used for purposes other than the production of permanent single crowns, inlays, onlays and veneers. Any deviation from these instructions for use can have negative effects on the chemical and physical quality of plastics made from VarseoSmile Crown<sup>plus</sup>.

### 3. Safety instructions

This product is produced and tested according to the most stringent quality standards. In order to ensure optimum further processing, please read the information contained in the instructions for use carefully.

The safety and precautions included in these instructions for use and safety data sheet shall apply to the handling of liquid resin and printed objects that have not been post-cured (objects in the “green condition”).

#### 4. Side effects and precautions

##### Precautions/Protection

It is essential that protective clothing be worn when handling this product. Safety goggles and nitrile gloves must be used. Further information on handling the product can be found in the safety data sheet and also downloaded from the BEGO Download Centre at [www.bego.com](http://www.bego.com). We cannot completely rule out adverse reactions (e. g. intolerance or allergies) to specific material components for all individuals. In such isolated cases, the user should discontinue use of the material.

**WARNING**

##### Contains:

Esterification products of 4,4'-isopropylidiphenol, ethoxylated and 2-methylprop-2enoic acid.  
Silanized dental glass, methyl benzoylformate, diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide.  
Total content of inorganic fillers (particle size 0.7 µm) is 30–50% by mass.

Hazard statements as per MSDS

- Causes skin irritation.
- May cause an allergic skin reaction.
- Causes serious eye irritation.
- May cause respiratory irritation.
- May cause long harmful effects to aquatic life.

Precautionary statements as per MSDS

- Avoid breathing mist/vapours/spray.
- Avoid release to the environment.
- Wear protective gloves/protective clothing/eye protection/face protection.
- Call a POISON CENTER/doctor if you feel unwell.
- If skin irritation or rash occurs: Get medical advice/attention.
- If eye irritation persists: Get medical advice/attention.
- Dispose of contents/container as per local and national regulations.

**Notice:** Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## 5. General information on handling

### Delivery

VarseoSmile Crown<sup>plus</sup> is supplied in seven colours according to the VITA\* classical shade system, in lightproof and sealed bottles.

Filling quantity:

- REF 41107 = 500 g, A1 Dentin
- REF 41117 = 250 g, A1 Dentin
- REF 41108 = 500 g, A2 Dentin
- REF 41118 = 250 g, A2 Dentin
- REF 41109 = 500 g, A3 Dentin
- REF 41119 = 250 g, A3 Dentin
- REF 41110 = 500 g, B1 Dentin
- REF 41120 = 250 g, B1 Dentin
- REF 41111 = 500 g, B3 Dentin
- REF 41121 = 250 g, B3 Dentin
- REF 41112 = 500 g, C2 Dentin
- REF 41122 = 250 g, C2 Dentin
- REF 41113 = 500 g, D3 Dentin
- REF 41123 = 250 g, D3 Dentin

### Storage

This product must be stored in the original sealed bottle, or in the cartridge at room temperature (approx. 22 °C) in a dark, dry place. It must be ensured that the temperature does not drop below +4 °C and does not exceed +28 °C! The minimum shelf life date printed on the product must be observed. Perfect processing cannot be guaranteed if materials which have exceeded their minimum shelf life date are used.

**Precaution:** Expected results cannot be guaranteed if materials which have exceeded their minimum shelf life date are used or if storage instructions are not followed.

The completely cured print objects must be stored at room temperature and protected from sources of light.

\* This symbol is a commercial designation/registered trademark of a company that is not part of the BEGO company group.

## 6. Processing requirements

### 1. Design

- Create the object (STL-file) using a commercial CAD software, which is intended for dental applications.
- When designing, observe the requirements for minimum wall thicknesses for finished restorations.

#### Notice:

##### Single crowns, inlays, onlays and veneers

Minimum wall thicknesses anterior teeth	1.0 mm
Minimum wall thicknesses posterior teeth	1.0 mm
Minimal wall thickness, cervical	1.0 mm
Minimal wall thickness, veneer shell	0.5 mm

### 2. Compatible 3D printers for additive manufacturing

#### Compatible 3D Printers

3D Printer Model	Printer Firmware	Nesting Software	Provider
Varseo	1.14 or higher	BEGO CAMCreator Print Version 1.14 or higher	BEGO
Varseo L	1.02 or higher		
Varseo S	1.14 or higher		
Varseo XS	2.6.8.24 or higher		
ASIGA MAX UV	2018-09-03	ASIGA Composer	ASIGA
SprintRay Pro 55	6.32.5	RayWare	SprintRay
SprintRay Pro 95	2.5.1	RayWare	SprintRay

### 3. Nesting & preparation for printing

- Import STL file
- Manual/Automatic rotation and placement
- Optimal orientation: horizontal, occlusal plane facing the build platform
- Manual/Automatic generation of supports

### 4. Necessary tools, equipment and materials for post-processing

- Stainless steel spatula
- Unheated ultrasonic bath
- Ethanol solution 96 %
- Spray bottle with 96 % ethanol solution
- Cutting wheel or side cutters (for support structure removal)
- Sandblaster 1.5 bar
- Glass bead blasting material 50 µm (e. g. Perlablast® micro, REF 46092/54302)

### 7. Processing

VarseoSmile Crown<sup>plus</sup> is compatible and usable with 3D printing system components from various providers. You can find these compatible system components on our website <https://www.bego.com/3d-printing/compatibility-overview/>

The printer settings can be found in the manual for use for the equipment.

**Precaution:** The material specifications have been validated using the software, printers, and process parameters specified in this document. Any unauthorized changes to the process, equipment, or software may result in a device that is out of specification. When using VarseoSmile Crown<sup>plus</sup>, users shall follow all instructions found in this document. The user shall also follow all manufacturer provided instructions for the use and maintenance of the associated processing devices (e. g. printers and post-cure devices).

The following instructions contain details of a validated workflow for the 3D printing process.

### Printing

Follow the instruction for use of the 3D printer manufacturer's to prepare the printing process. The ideal working temperature of VarseoSmile Crown<sup>plus</sup> is in the temperature range between 18 °C and 28 °C.

Before use, the resin must be homogeneous. Before the first use, the material has to be shaken well about 2 min. When decanting, make sure that the printing resin is exposed to daylight for as short a period of time as possible. Mix the resin in the cartridge/resin tank if a transparent layer is visible on the surface.

For further processing – selecting the resin, material parameter, setting up the print job – as part of the printing process, follow the respective printer instructions for use.

### Cleaning and preparation for post-curing

On completion of printing, the print objects are released from the build platform using the spatula. The print object should be cleaned in two steps with ethanol (96 %) using an ultrasonic bath.

**Precaution:** Never fill ethanol directly into the ultrasonic bath; place it in the recommended container (REF 19621) in the ultrasonic bath filled with water. Use an explosion-proof ultrasonic bath.

1. Clean the print object for **3 min** in a reusable ethanol solution (96 %) using an **unheated** ultrasonic bath.
2. The precleaned object must be cleaned thoroughly for **2 min** using a fresh ethanol (96 %) solution with the aid of an **unheated** ultrasonic bath. The print object is then removed from the ethanol bath and sprayed with additional ethanol (96 %) in order to fully rinse off any remaining resin residue.

**Tip:** Resin residues can also be removed using a brush soaked in ethanol (96 %).

**Precaution:** The entire cleaning process should not take longer than 5 min as this could otherwise have a detrimental effect on the printed objects (swelling of the object with ethanol).

After cleaning, the print object is dried using compressed air under an extraction unit. If there is liquid resin still adhering to the surface of the object, this can be completely removed by spraying again with ethanol (96 %) and re-drying.

### Preparation for post-curing

- Remove the support structure with the help of a cutting wheel or side cutters.
- Remove the white layer using a glass bead blasting material 50 µm (e. g. Perlablast® micro, REF 46092/54302) at a maximum blasting pressure of 1.5 bar.
- Check for fit and finish the objects completely. Finishing and countouring can be performed using carbide cutter or diamond grinding stones.

### Post-curing process

The final properties of the printed object depend on the post-curing process. Please note the assignment of the light curing device to the 3D printer of the approved system components.

The post-curing of the object is done without use of a model, then allow to cool for 3–5 minutes or until the object feels cool.

Post-curing			
3D Printer	Exposure Cycles	Light-curing Device	Additional Information
BEGO Varseo	BEGO-Otoflash	2 x 1,500 flashes	Turn object between the exposure cycles
BEGO Varseo L			
BEGO Varseo S			
BEGO Varseo XS	HiLite-Power*	2 x 90 seconds	
ASIGA MAX UV			
SprintRay Pro 55	SprintRay ProCure	2 x 20 min at 20 °C	
SprintRay Pro 95	SprintRay ProCure	2 x 20 min at 20 °C	

**Note:** When employing the BEGO Otoflash, use the protective gas function. This results in a further reduction of the already low remaining monomer content. To do so, set the protective gas function to switch position 1. Details can be found in the instruction manual for the post-curing device.

The times given only apply to regularly maintained equipment that guarantees a corresponding light intensity.

\* This symbol is a commercial designation/registered trademark of a company that is not part of the BEGO company group.

**Precaution:** If there is an interruption or failure in the post-curing device cycle, the printed object should be not used until it has cured under a full cycle. Check the post-curing device manual for use for how to properly resolve the post-curing device condition and then repeat the post-curing cycle with the printed objects.

### Supplementing of printed objects

Defects (e. g. missing contact points) can be supplemented with the resin or with commercially available composite veneering materials.

### Supplementing with resin

- Blast the areas to be supplemented with aluminum oxide – pressure 1.5 bar and particle size 110 µm.
- Put some liquid VarseoSmile Crown<sup>plus</sup> on the object.
- Put it under light for a short time (5 flashes) in the BEGO Otofash
- If more material has to be applied, put additional liquid VarseoSmile Crown<sup>plus</sup> on the last layer and then put it again under light for a short time (5 flashes) in the BEGO Otofash.
- Final curing of the Objects by using the BEGO Otofash for 2 x 1500 flashes or HiLite Power for 2 x 90 s.

### Supplementing with veneering material

The printed objects can be supplemented with commercially available composite veneering materials (e. g. VITAVM LC\*, Vita Zahnfabrik, und VITAVM LC flow\*, Vita Zahnfabrik). The manufacturer's instructions for use must be observed.

### Polish

Polish the surface of the objects with pumice stone and polishing compound. Avoid overheating of the restorations during polishing. Optimal surface quality is achieved by polishing after post-curing.

**Tip:** Optionally, the surface of the objects can be coated with light-curing glaze (e. g. Vita ENAMIC GLAZE \*, Vita Zahnfabrik or GC OPTIGLAZE \*, GC). Pay attention to the manufacturer's instructions for use.

### Individualization (optional)

Instead of polishing the surface, the post-cured objects can also be individualized with composite stains or veneering materials.

\* This symbol is a commercial designation/registered trademark of a company that is not part of the BEGO company group.

### **Individualization with composite stain**

Commercially available composite stain materials (e. g. VITA ENAMIC STAIN\*, Vita Zahnfabrik and GC OPTIGLAZE Color\*, GC) can be used for individualization of the post-cured VarseoSmile Crown<sup>plus</sup> restoration. The manufacturer's instructions for use must be observed.

### **Individualization with veneering materials**

Commercially available composite veneering materials (e. g. VITAVM LC\*, Vita Zahnfabrik and VITAVM LC flow\*, Vita Zahnfabrik) can be used for individualization of the restoration. The manufacturer's instructions for use must be observed.

### **Veneer shells**

Commercially available light-curing veneering composites can be used to attach the post-cured veneer shells to the metal framework (e. g. VITAVM LC, VITAVM LC flow, Vita Zahnfabrik).

### **Basic work steps**

- Condition the surface of the metal framework.
- Cover the metal framework with opaquer.
- Sandblast the inside of the veneer shell and clean the veneer shell from dust.
- Inner surface conditioning of the veneer shell. Observe the information of the veneering material manufacturer in section "Individualization of composite frameworks" or "Individualization of artificial teeth".
- Applying the veneering material to the inside of the veneer shell or the metal framework.
- Positioning of the veneer shell on the metal framework and remove the excess with a brush or an instrument.
- Carry out the polymerization according to the instructions of the veneering material manufacturer.
- Finally, the surface of the veneer shell should be polished or customized with composite stains (see section "**Polishing**" or "**Individualization**").

\* This symbol is a commercial designation/registered trademark of a company that is not part of the BEGO company group.

## 8. Cleaning in the dental laboratory and dental practice

Fully cured objects made from VarseoSmile Crown<sup>plus</sup> can be easily cleaned and disinfected. Steam cleaning (e. g., with Triton SLA) is possible. Disinfection in the immersion bath (e. g. ethanol 96 % or MD 520\* impression disinfectant, Dürer Dental Co.) is also possible. Follow manufacturer's instructions.

## 9. Luting

### Attachment to adhesive Titan abutment/Titan abutment

Commercially available adhesive luting materials (e. g. Panavia V5\*, Kuraray Noritake) can be used for fastening the Crown to the adhesive Titan abutment / Titan abutment. The manufacturer's instructions for use must be observed.

### Observe the following points when fixing the crowns on titanium adhesive bases:

- After the restoration is attached to the implant, a Teflon band is pressed into the screw channel to close the screw channel.
- The screw channel is then closed with a composite (e. g. Tertric EvoCeram Bulk Fill\*, Ivoclar Vivadent).
- The objects can be polished with commercially available composite polishers.

### Attachment to tooth stump

The finished permanent restorations can be attached using self-adhesive cements (e. g. RelyX Unicem\*, 3M Espe) or composite cement with a primer (e. g. Variolink Esthetic DC\* and Monobond Plus\*, Ivoclar Vivadent). Observe the instructions for use of the luting agent.

## 10. Disposal

The cured, separated material (base plate, support structure) can no longer be used. Cured material can be disposed of as domestic waste. Unused resin or ethanol used for cleaning with resin residues must be disposed of via the local waste disposal authority or a hazardous waste collection point stating the safety data sheet.


\* This symbol is a commercial designation/registered trademark of a company that is not part of the BEGO company group.

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**BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG**