dima Print Denture Base Try-in

dima Print Denture Base Try-in

- DE** Gebrauchsanweisung
- GB** Instructions for use
- FR** Mode d'emploi
- ES** Instrucciones de uso
- IT** Istruzioni per l'uso
- PT** Instruções de uso
- NL** Gebruiksaanwijzing

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Stand: 2020-10

GB Instructions for use dima Print Denture Base Try-in

Indications for Use

dima Print Denture Base Try-in is a light-curable resin indicated for fabrication of Try-in denture before deciding to make a final denture. The material is an alternative to traditional heat-curable and auto polymerizing resins to make a Try-in denture.
The fabrication of Try-in denture with dima Print Denture Base Try-in requires following components: digital denture files based on a digital scan of an impression, a stereolithographic additive printer, and curing light equipment. The use of Try-in denture is to check out the esthetic and functional aspects at the chair side for a short amount of time.

Requirements

- Digitale model file with minimum thickness 2 mm; STL format file
- Additive printer and its operation software;

	Additive Printer	Operation Software	Provider
	cara Print	cara Print CAM	Kulzer

Curing Equipment	Provider
Hilite power 3D	Kulzer

Specific Manufacturing Considerations

1. Digital denture model file

- File format: STL file
- Minimum thickness of denture model: ≥ 2 mm
- File size should be uploadable in the 3D printer operation software.

2. Additive printer

- Hardware
 - Wavelength: 385 nm or 405 nm
 - Light source
 - Stereolithographic (SLA) method; laser with $25\text{ mW} < X < 250\text{ mW}$
 - Digital Light Processing (DLP) method; high power LED or lasers
 - Build Volume: $> 103 \times 58 \times 130$ mm (Least fit one arch)
 - Build Path: line drawing path or surface layer drawing path
- Features of Operation Software
 - STL file import
 - Automatic rotation and placement
 - Layer slicer for path inspection
 - Auto and manual generation of supports

2.3 Printing Parameter

Printer Model	Layer Thickness (micron)	Recommended orientation angle (degree)	Support point size (mm)	Support density (mm)
cara Print	50, 100	20–40	1–2	min. 1

2.4 Environmental Conditions

- Temperature: $18\text{--}30^{\circ}\text{C}$ ($64\text{--}86^{\circ}\text{F}$)
- Relative Humidity: 30–90%

2.5 Cleaning Kit

Rinse bath and tubs, flush cutter, paper towel, squeeze bottle for isopropyl alcohol, Scraper

3. Recommended Curing light equipment (Post curing units)

3.1 Flood Type Curing Equipment

Manufacturer / Model	Curing Chamber	Supply voltage (V / Hz)	Lamp power (W)	Light intensity (mW/cm ²)	Lamp wavelength (nm)	Curing time (min)
Kulzer / Hilite power 3D	Required	100,115,230 / 50–60	200	N/A	390–540	20

3.2 Accessories

- USP Grade glycerin
- Transparent glass container and 2 glass plates
- Heat-protective gloves and silicone coated stainless steel tong
- Thermocouple

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Manufacturer:
DENTCA, INC.
357 VAN NESS WAY, STE 250
TORRANCE, CA 90501
USA

Importer:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)
+49 800 4372522

Distributed in USA /
Canada exclusively by:
Kulzer, LLC
4315 South Lafayette Blvd.
South Bend, IN 46814-2517
1-800-431-1785

Caution: Federal law restricts
this device to sale by or on the
order of a dental professional.

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EU Representative:
MT Promed Consulting GmbH
Altenhofstr. 80
66386 St. Ingbert
Germany

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Stand: 2020-10

4. Notification

The device specifications have been validated using the software, printers, and process parameters specified in this document. Any other printers, operation software and post-printing processes will be outside of the device specifications and the FDA clearance. Users shall follow this document to use the device.

Warnings:

- dima Print Denture Base Try-in resin contains polymerizable monomers which may cause skin irritation (allergic contact dermatitis) or other allergic reactions in susceptible persons. If contact with skin, wash thoroughly with soap and water. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical assistance.
- Avoid inhalation or ingestion. High vapor concentration can cause headache, irritation of eyes or respiratory system. Direct contact with eyes may cause possible corneal damage. Long-term excessive exposure to the material may cause more serious health effects. Monitor air quality per OSHA standards.
- Eye Contact: Immediately flush eyes with plenty of clean water for at least 20 minutes, and consult a physician. Wash the contacted area thoroughly with soap and water.
- Inhalation: In case of exposure to a high concentration of vapor or mist, remove person to fresh air. Give oxygen or artificial respiration as required.
- Ingestion: Contact your regional poison control center immediately.
- BURN HAZARD: GLYCEROL BATH CAN REACH TEMPERATURES OF 90°C ($\sim 200^{\circ}\text{F}$) AND LEAD TO SEVERE BURNS. Only trained users should perform the glycerol curing step with caution and appropriate PPE. We also recommend placing a warning label on the window of the cure unit to alert all lab users to the potential hazard.

Precautions:

- When washing the printed Try-in denture with solvent or grinding the Try-in denture, it should be done in a properly ventilated environment with proper protective masks and gloves.
- Store dima Print Denture Base Try-in resin at or below $15\text{--}25^{\circ}\text{C}$ ($60\text{--}77^{\circ}\text{F}$) and avoid direct sunlight. Keep container closed when it is not in use. Product shall not be used after expiration date.
- Expired or unused dima Print Denture Base Try-in should be completely cured or polymerized prior disposal.

Adverse Reactions:

- Direct contact with the uncured resin may induce skin sensitization in susceptible individuals.
- Proper ventilation and personal protective equipment should be used when grinding printed Try-in denture as the particulate generated during grinding may cause respiratory, skin and eye irritation.

Procedure to Fabricate the Try-in Denture

- Printing Preparation
 - Open the 3D printer cover and fill the resin tank or tray of the printer with dima Print Denture Base Try-in resin up to the required filling line by manufacturer.
(When filling the resin into the resin tank or tray, gloves and mask should be used.)
 - Close the printer cover.
- Printing
 - Load the denture model file to be printed in printer operation software which printer manufacturer recommended.
 - Use auto-orientation or manual orientation to find its optimal position for printing. The recommended orientation by printer provider is a tilted orientation such as space diagonal from 20 to 90 angle. If auto-orientation is not satisfied, rotate to make optimal position.
 - Generate support sticks on the denture model. The support is not enough, add supports on the denture model.
 - Use layout tools of the software to move the denture model within the imaginary build platform to prevent the overlapping between models.
 - Start printing.
- Cleaning
 - Detach the printed Try-in denture from the build platform.
 - Use a small flush cutter to remove the support sticks from the printed Try-in denture.
 - Wash the Try-in denture with isopropyl alcohol.
 - Use air blowing to dry the Try-in denture or dry it at room temperature under ventilation system or open area.
- Post Processing of Try-in denture
 - Add dima Print Denture Base Try-in material manually on the lingual ridge area to thicken and cure it. It is recommended to make it thicker than 4.0mm for maxillary and 4.5mm for mandibular.
 - Cure the printed Try-in denture by sinking into glycerin container for the required curing time under recommended post-curing unit. For half of the post-curing time tissue side up and for another half of the time tissue side down.
Note: It is recommended the glycer

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