Instructions for Use

DENTURE TEETH RESIN

Prepared: 09/24/2020 PRNT-107234 Rev 01

The following Instructions For Use are for Formlabs biocompatible photopolymer Denture Teeth Resin. Basic information about safety and environmental concerns are also included. For more detailed safety and environment information please refer to the Safety Data Sheet, available at dental.formlabs.com. For further information regarding the use of the material, please contact Formlabs.

1. Introduction And Indications For Use

INDICATIONS FOR USE

Denture Teeth Resin is a light-curable polymerizable resin to fabricate, by additive manufacturing, preformed denture teeth to be used in a denture. The fabricated tooth is an alternative to preformed plastic tooth for denture.

The fabrication of denture tooth with Denture Teeth Resin requires digital denture tooth files instead of physical moulds, a Formlabs SLA 3D printer, and light curing equipment.

2. Specific Manufacturing Considerations

NOTIFICATION

The device specifications have been validated using the software, printers, and process parameters specified in this document. Any other printers, operation software and post-printing processes will be outside of the device specifications and the FDA clearance. Users shall follow this document to use the device.

REQUIREMENTS

Use dedicated accessories for Denture Teeth Resin. For biocompatibility, Denture Teeth Resin requires a dedicated resin tank, build platform, Form Wash, and finishing kit, which should not be mixed with any other resins.

RECOMMENDED 3D PRINTER AND PRINTING PARAMETERS

- a. Hardware: Formlabs SLA 3D printer
 - Laser wavelength: 405 nm
 - Z axis layer thickness: 50 μ m
 - Laser power: 25 mW < X < 250 mW
 - Build volume: > 70 x 50 x 150 cm (at least fit one arch)
 - Laser spot size (XY resolution): < 160 μ m
 - Build speed: 2 - 3 cm/hr at 50 μ m
- b. Software: Formlabs PreForm
 - STL file import
 - Manual/Automatic rotation and placement
 - Layer slicer for path inspection
 - Auto-mesh repair
 - Manual/Automatic generation of supports
- c. Printing Parameters
 - Layer thickness: 50 μ m
 - Optimal orientation: 20° - 40° tilted orientation
 - Support point size: 0.4 - 1.0 mm
 - Support density: 0.7 - 1.5
 - Minimum thickness of 2 mm
- d. Recommended Post-Processing Equipment:
 - Formlabs Form Wash
 - Isopropyl alcohol (IPA) $\geq$ 90%
 - Formlabs Form Cure
 - USP grade glycerol ($\geq$ 99.5% purity, CAS# 56-81-5)
 - Transparent glass container (recommend Pyrex* Basics 23oz dish: 6" x 6" x 2.5")

*Pyrex is a registered trademark not affiliated with Formlabs, Inc.

3. Hazards And Precautions

WARNINGS

1. Denture Teeth Resin contains polymerizable monomers which may cause skin irritation (allergic contact dermatitis) or other allergic reactions in susceptible persons. If resin contacts skin, wash thoroughly with soap and water. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical assistance.
2. Avoid inhalation or ingestion. High vapor concentration can cause headache, irritation of eyes or respiratory system. Direct contact with eyes may cause possible corneal damage. Long-term excessive exposure to the material may cause more serious health effects. Monitor air quality per OSHA standards.
3. **BURN HAZARD:** GLYCEROL BATH CAN REACH TEMPERATURES OF 90°C (~200°F) AND LEAD TO SEVERE BURNS. Only trained users should perform the glycerol curing step with caution and appropriate PPE. We also recommend placing a warning label on the window of the cure unit to alert all lab users to the potential hazard.
4. **Eye contact:** Immediately flush eyes with plenty of clear water for at least 20 minutes, and consult a physician. Wash the contacted area thoroughly with soap and water.
5. **Inhalation:** In case of exposure to a high concentration of vapor or mist, remove the person to fresh air. Give oxygen or artificial respiration as required.
6. **Ingestion:** Contact your regional Poison Control Center.

PRECAUTIONS

1. When washing the printed part with solvent or grinding the denture teeth, it should be in a properly ventilated environment with proper protective masks and gloves.
2. Expired or unused Denture Teeth Resin shall be disposed in accordance with local regulations.
3. Glycerol and isopropanol shall be disposed of in accordance with local regulations.

ADVERSE REACTIONS

1. Direct contact with the uncured resin may induce skin sensitization in susceptible individuals.
2. Proper ventilation and personal protective equipment should be used when grinding denture teeth as the particulate generated during grinding may cause respiratory, skin and eye irritation.

4. Procedure To Fabricate Dentures

A. PRINTING AND POST-PROCESSING

1. Select the Denture Teeth Resin shade based on prescription. **Note:** Use a dedicated resin tank for different shades.
2. **Shake cartridge:** Prior to use, shake cartridge well. Color deviations and print failures may occur if the cartridge is shaken insufficiently.
3. **Set up:** Insert resin cartridge into a compatible Formlabs 3D printer.
4. **Printing:**
 - a. Prepare a print job using PreForm software. Import the denture teeth STL file. Orient and generate supports. For recommendations on print orientation and support placement, see the detailed application guide at dental.formlabs.com.
 - b. Send print job to printer. Begin print by selecting a print job from the print menu. Follow any prompts or dialogs shown on the printer screen. Printer will automatically complete print.
5. **Part removal:**
 - a. Remove the build platform from the printer.
 - b. Printed parts can be removed from the build platform before or after cleaning in a Form Wash. To remove, wedge the part removal tool under the printed part, and rotate the tool. For detailed techniques visit support.formlabs.com.
6. **Rinsing:** Wash the printed parts in IPA using a Formlabs Form Wash for 10-20 minutes.
7. **Drying:**
 - a. Remove parts from IPA and either leave to air dry at room temperature, or dry using a source of compressed air.
 - b. Inspect parts to ensure that the denture teeth are fully washed and dried with no particles or residues before proceeding to subsequent steps.
8. **Support removal:**
 - a. Remove supports using a cutting disk and handpiece, or with other part removal tools.
 - b. Carefully smoothen the supported surfaces of the denture teeth using a bur and handpiece.

9. Teeth Post-Curing:

Attention: For printed dentures, do not post-cure printed denture teeth prior to bonding with denture base.

- a. For conventional fabrication of preformed teeth, fill a glass container with glycerol, and place it on the turntable inside a Form Cure. 6" diameter Pyrex* containers are recommended.
- b. Set the Form Cure to 80°C and wait for 2 hours until the glycerol fully warms up.
- c. Place the teeth in the glycerol, ensuring that it is fully submerged.
- d. Cure for 30 minutes at 80°C.
- e. Flip the teeth and cure for an additional 30 minutes at 80°C.
- f. Use the preformed teeth to fabricate a denture by conventional denture process.

B. DENTURE FABRICATION USING PRINTED DENTURE BASE AND TEETH

1. Bonding denture teeth to denture base:

- a. Before post-curing, remove any material from the teeth and tooth sockets on the denture base, which may impede tooth placement.
- b. Place the printed teeth with Denture Teeth Resin into the corresponding tooth sockets on the printed denture base and check teeth fitting.
- c. Remove the teeth to prepare for bonding.
- d. Apply a small amount of Denture Base Resin into the tooth sockets using a brush, press the teeth into the corresponding socket and hold firmly in place.
- e. Tack the teeth into place using a handheld UV-light source, applying light to both the lingual and buccal/facial surfaces.
- f. If necessary, apply a small amount of Denture Base Resin to the interproximal surfaces or elsewhere to ensure a smooth, complete bond between the teeth and base.
- g. Repeat above steps until all teeth are bonded to the denture base.
- h. Check the occlusion of the assembled denture against the opposing denture/model before post-curing.
- i. Add Denture Base Resin manually on the lingual ridge area to thicken and cure it. It is recommended to make it thicker than 4.0mm for maxillary and 4.5mm for mandibular.

*Pyrex is a registered trademark not affiliated with Formlabs, Inc.

2. **Post-curing:**

Attention: Only post-cure the denture base after fully assembling with the denture teeth.

- a. Fill a glass container with glycerol, and place on the turntable inside a Form Cure. 6" diameter Pyrex* containers are recommended.
- b. Set the Form Cure temperature to 80°C and wait for 2 hours until glycerol fully warms up.
- c. Place the assembled denture in the glycerol, ensuring that it is fully submerged.
- d. Cure for 30 minutes at 80°C.
- e. Flip the assembled denture and cure for an additional 30 minutes at 80°C.

3. **Finishing:**

- a. Polish the final denture using conventional methods as for traditional dentures. It is recommended to wet polish using a rag wheel and pumice, using low speeds initially to smoothen the surfaces without removing characterization, followed by a quick polish at high speeds using acrylic high-shines.

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DENTCA, Inc.