## 9880906 - Precision Investment for Model Casting Technique

EN

### INSTRUCTIONS FOR USE

Please ensure that the equipment is free from any residue including cleaning agents, which might have a negative effect on the casting result. Use a separate spatula and mixing bowl when working with investments, which should be rinsed after each use. Keep all equipment used for investments separate to items used for stone and plasterwork. Always use a fine filtration dust mask when weighing powder or devesting the mould.

**Application:** Chrome investment suitable for gel and silicone duplication in the traditional heating process, silicone duplication by rapid burnout and light cure wax technique in the traditional heating process.

**Technical data:** Working time: 4 - 6 minutes; Working temperature: 20 - 23°C powder and liquid (ideally at 21°C in a temperature chamber, in case of gel duplication 24 - 25°C).

**Mixing ratio:**

|                              | Powder | Total liquid | Volume of Mixing Liquid for Precision Investment | Volume of demineralized water |
|------------------------------|--------|--------------|--------------------------------------------------|-------------------------------|
| gel duplication              | 100 g  | 19 ml        | 13.3 ml (70%)                                    | 5.7 ml (30%)                  |
| silicone duplication         | 100 g  | 21 ml        | 16.8 ml (80%)                                    | 4.2 ml (20%)                  |
| imbedding                    | 100 g  | 21 ml        | 10.5 ml (50%)                                    | 10.5 ml (50%)                 |
|                              | 400 g  | 84 ml        | 42 ml (50%)                                      | 42 ml (50%)                   |
| light cure wax (see point 9) | 100 g  | 21 ml        | 16.8 ml (80%)                                    | 4.2 ml (20%)                  |
|                              | 600 g  | 126 ml       | 101 ml (80%)                                     | 25 ml (20%)                   |

Room temperature, air humidity and mixer adjustment can influence the results.

**Recommendations and tips regarding expansion:** Changes to volume of liquid can influence expansion. More mixing liquid – greater expansion, more distilled water - less expansion. A change to the total liquid used (up to 4ml) will result in either a thicker mix – bigger expansions (wider casting) thinner mix – low expansions (tighter casting).

**Processing:** Recommendation: A silicone mould with a shore hardness of 17 to 22 should be used without stabilization or adhesive sleeve thus enabling the investment to expand freely.

**Model preparation using silicone duplication:** Put the powder into the mixing bowl and weigh (100 g = 21 ml). Add the (mixed up) total liquid (begin timing 20 minutes). Mix vigorously by hand for 15 seconds. Mix under vacuum for 45 seconds; mixing speed 250 Rpm. Fill in the investment evenly from one side only at low vibration level. The thinnest part of the model should be at least 1 cm thick. Demould the model after 20 minutes and dry at 140°C for 20 minutes. Do the wax up. Fix sprues and cone. Mix the investment for imbedding, start timing. Pour the investment in without vibrating. After 20 minutes place the mould into the furnace for rapid burn out.

**Model preparation using gel duplication (for example with duplicating gel from Henry Schein):** The gel mould should be dry and lukewarm before pouring into the investment. Put the powder into the mixing bowl and weigh (100 g = 19 ml). Add the (mixed up) total liquid (Begin timing 30 minutes). Mix vigorously by hand for 15 seconds. Mix under vacuum for 60 seconds; mixing speed 250 rpm. Fill in the investment evenly from one side at low vibration level. The thinnest part of the model should be at least 1 cm thick. Demould the model after 30 minutes and then dry it at 180°C for 30 minutes. Dip for 2 seconds into the dipping hardener. Do the wax up. Fix sprues and cone.

**Overbedding:** Place the model onto the hole on the bottom of the muffle. Fix the model at the bottom of the muffle using adhesive wax (higher melting point). Place the thinner side of the mould ring on the bottom of the mould. Mix the investment for imbedding (Begin timing 20 minutes). Pour the investment into the mould without vibrating. After 20 minutes put the mould into the furnace for rapid burn out.

**Heating / Preheating:** Put the mould into the preheated furnace onto a punched or a coarse milled ceramic plate with the cone facing the bottom. In case of gel duplication we recommend using the traditional heating process.

**Traditional heating process:** Place the mould into the cold furnace, as soon as possible after 20 minutes mixing process. Heating rate: up to 20°C/min, (holding times are not necessary, except for the light curing technique). End temperature (according to the alloy used) hold for at least 45 minutes. We recommend an end temperature of 850°C.

**Rapid burn out:** Place the mould into the hot furnace 20 minutes after the beginning of the mixing process – at a maximum temperature of 850°C for at least 45 minutes. If necessary, you can heat up until end temperature.

**Casting:** After a holding time of at least 45 minutes at end temperature the casting can begin in accordance with the alloy manufacturer's instructions. If several moulds are going to be heated up in the furnace, the holding time has to be extended by 10 minutes per mould.

**Cooling down:** Cool down mould slowly to room temperature.

**Special considerations of the light curing technique:** Humidify the light curing wax surface with the surface tension release agent before pouring in the investment. Imbedding (see table). Traditional heating process. Holding time: 30 minutes at 350°C. End temperature: 850°C (hold for at least 45 minutes).

**Health warning:** Investments contain quartz! Do not breathe dust. Danger of lung diseases (silicosis or cancer).

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**HENRY SCHEIN U.K. HOLDINGS LTD.**