


# AcryPlus

*We take care of your smile*

DENTS ARTIFICIELLES EN RÉSINE - DIENTES ARTIFICIALES EN RESINA



**Échelle couleurs - Escala de colores**

Acry Plus	Acry Plus V
1A - 2A - 3A - 4A	A1 - A2 - A3 - A3,5 - A4
2B - 4B - 5B - 6B	B1 - B2 - B3 - B4
1C - 2C - 3C - 4C - 6C	C1 - C2 - C3 - C4
1D - 4D - 6D	D2 - D3 - D4
1E - 2E - 3E	

**T1**

(1) 39,4  
(2) 9,2  
(3) 11,5  
(4) 7,2

**T2**

(1) 42,5  
(2) 9,5  
(3) 11,3  
(4) 7,7

**T3**

(1) 41,9  
(2) 10,4  
(3) 12,4  
(4) 8,2

**T4**

(1) 43,0  
(2) 10,4  
(3) 12,4  
(4) 8,2

**T5**

(1) 46,4  
(2) 10,5  
(3) 13,0  
(4) 8,4

**T6**

(1) 46,5  
(2) 11,0  
(3) 13,0  
(4) 8,5

**T7**

(1) 48,2  
(2) 10,7  
(3) 13,4  
(4) 8,5

**T8**

(1) 43,8  
(2) 9,6  
(3) 11,4  
(4) 7,6

**T9**

(1) 41,4  
(2) 9,0  
(3) 11,0  
(4) 7,6

**T16**

(1) 45,5  
(2) 10,2  
(3) 11,5  
(4) 8,3

**T17**

(1) 46,1  
(2) 10,2  
(3) 11,5  
(4) 8,3

**T18**

(1) 46,7  
(2) 10,3  
(3) 11,5  
(4) 8,5

**T19**

(1) 48,1  
(2) 10,7  
(3) 12,4  
(4) 8,9

**P1**

(1) 41,0  
(2) 8,2  
(3) 9,7  
(4) 7,8

**P2**

(1) 44,5  
(2) 8,6  
(3) 10,2  
(4) 8,2

**P3**

(1) 43,0  
(2) 9,4  
(3) 10,7  
(4) 8,2

**P4**

(1) 45,2  
(2) 10,6  
(3) 11,7  
(4) 8,7

**P5**

(1) 48,3  
(2) 10,5  
(3) 13,0  
(4) 8,8

**P7**

(1) 51,0  
(2) 12,0  
(3) 14,0  
(4) 9,4

**P8**



(1) 45,5  
(2) 12,0  
(3) 14,0  
(4) 9,4

**P9**



(1) 49,4  
(2) 9,3  
(3) 11,0  
(4) 9,0

Code produit Código producto	
0120010	ACRY PLUS 6 UPPER
0120011	ACRY PLUS 6 LOWER
0120012	ACRY PLUS 8 UPPER
0120013	ACRY PLUS 8 LOWER
0120030	ACRY PLUS N.M. 8 UPPER
0120031	ACRY PLUS N.M. 8 LOWER
0120014	ACRY PLUS SET 28
0120015	LIVE MOULD CHART ACRY PLUS
0120016	SHADE GUIDE ACRY PLUS

0120020	ACRY PLUS V. 6 UPPER
0120021	ACRY PLUS V. 6 LOWER
0120022	ACRY PLUS V. 8 UPPER
0120023	ACRY PLUS V. 8 LOWER
0120032	ACRY PLUS V. N.M. 8 UPPER
0120033	ACRY PLUS V. N.M. 8 LOWER
0120024	ACRY PLUS V. SET 28
0120025	LIVE MOULD CHART ACRY PLUS V.
0120026	SHADE GUIDE ACRY PLUS V.




Disponibles dans les teintes Acry Plus V  
Disponibles en los colores Acry Plus V



Acry Plus comply with ISO 22112:2005 - Type 1 - Type 2

**Tableau d'articulation**  
**Tabla de articulación**



**G:** Géométrique - *Geométrico*  
**NM:** Neuro musculaire - *Neuromuscular*

P 1	X 1 - X 3	Z 3	
P 2	X 2 - X 5	Z 3	Z 14 - Z 17
P 3	X 3 - X 5	Z 3 - Z 5	Z 14
P 4	X 5 - X 6	Z 5	Z 14
P 5	X 7 - X 8	Z 4	Z 20 - Z 30
P 7	X 8 - X 9	Z 6	Z 20 - Z 30
P 8	X 5	Z 5	Z 17
P 9	X 7	Z 4 - Z 5	Z 20 - Z 30
T 1	X 3	Z 3	
T 2	X 3	Z 3 - Z 5	Z 14
T 3	X 2	Z 2	Z 14
T 4	X 4	Z 2	Z 14
T 5	X 10	Z 4	Z 17
T 6	X 7	Z 4	Z 17
T 7	X 8	Z 4	Z 20 - Z 30
T 8	X 5	Z 3 - Z 5	Z 14
T 9	X 3	Z 3 - Z 5	Z 14
T 16	X 5 - X 6	Z 5	Z 17
T 17	X 5 - X 7	Z 5	Z 17
T 18	X 5 - X 7	Z 4	Z 17
T 19	X 8 - X 10	Z 6	Z 20 - Z 30
V 1	X 3	Z 3 - Z 5	
V 2	X 5	Z 5	Z 17
V 3	X 7	Z 5	Z 17
V 4	X 1 - X 3	Z 3	Z 14
V 5	X 5	Z 5	Z 14 - Z 17
V 8	X 7 - X 10	Z 4	Z 20 - Z 30
V 9	X 8 - X 9	Z 6	Z 20 - Z 30

**Z2** (1) 29,9 (2) 31,5 (3) 9,0 (4) 8,5




**Z3** (1) 29,4 (2) 31,3 (3) 8,9 (4) 8,2



**Z5** (1) 30,5 (2) 33,2 (3) 9,1 (4) 8,8



**Z4** (1) 31,0 (2) 33,8 (3) 9,5 (4) 9,0



**Z6** (1) 33,6 (2) 36,3 (3) 9,8 (4) 9,1

